

## Specification Sheet - Washable Isolation Gown

- Universal (One Size)
- 75 washes at infection control temperatures
- Comfortable, quick dry stretch cuffs safely cover the wrist area
- Wrap-around design for enhanced coverage with tie fastenings
- Carbon yarns reduce static build-up
- Breathable fabric for comfort
- Designed with fewer seams for durability
- Stain-resistant properties
- Fluid Repellent fabric
- Widely used in Canada & USA

Description:	Washable Isolation Gown (Non-reversible)
Product Code:	GOWN/PPE
Fabric:	99% polyester, 1% Electrostatic Dissipative Carbon Yarns
Fabric weight:	105 gsm
Dimensions (cm):	Centre back length: 123cm, Sweep: 153cm, Cuff: 8.9cm
Colour:	Pale Yellow
Testing:	Fluid Repellent Fabric tested to AAMI Level 2 FDA registered
Laundering:	Wash at Thermal Disinfection Temperature Do not use bleach Do not exceed a maximum temperature of 90°C Tumble dry on a medium heat Remove load immediately once goods are dry 75 washes at recommended care instructions

*The washing process should have a disinfection cycle in which the temperature of the load is either maintained at 65°C for not less than ten minutes or 71°C for not less than three minutes when thermal disinfection is used – (Department of Health - Choice Framework for local Policy and Procedures 01-04 – Decontamination of linen for health and social care: Management and provision)*

**Please Note:** All measures and weights may vary by plus or minus 5% and are based on new, unwashed products

---