


S10-058-00100:
UPPER CONTROLLER



3D1200300000:
2A CHARGER

S01-046-00200:
BRACKETRY AS'Y


OPTIONAL EXTRAS:

- Z23014-08 EXTENDED TELESCOPIC BAR (FITS 20" TO 22" WHEELCHAIR SEAT WIDTH)
- Z10700-01 RHS PARKING BRAKE
- Z10750-01 LHS PARKING BRAKE
- C38-046-00100 : CARRY BAG




S18-58-00100:
WIRING LOOMS
(COMPLETE)

S01-046-00100:
TELESCOPIC BAR
(INC. 2
BARS &
KNOB)



S17-058-00200:
BATTERY PACK
(INCLUDES COVER,
BASE AND LEADS.
DOES NOT INCLUDE
BATTERIES)

S26-058-00100:
BATTERY TRAY AS'Y



C31-046-00200:
LIFT UP STRAP



3D1000100000:
PG S-DRIVE 45A



C30-001-10201:
M8 X 20 KNOB

C30-046-00300
QUICK REL. SCREW

C01-046-01500:
MOTOR HANDLE

MOTOR
(see note below)

MOTOR BRUSH
(see note below)

C09-005-02400:
MOTOR BRUSH CAP

C06-046-00100:
WHEEL C/W TYRE

MOTORS AND CONTROLLERS:

- C09-046-00200: 150W MOTOR
- C09-046-00600: 200W MOTOR (HEAVY DUTY)

MOTOR BRUSHES:

- C09-046-00900 : 6 x 8mm (FOR PWCPP009)
- C09-002-00101 : 6 x 10mm (FOR PWCPP009HD)



REV No.	DETAILS	DATE	DRAWN
A	ISSUE FOR PARTS	31/1/14	PMK
B	ISSUE FOR PARTS	16-8-15	DC
C	ISSUE FOR PARTS	06-07-16	PK

drive™
 Ainley's Industrial Estate, Eiland HX5 9JP
 T: +44 (0) 1422 314488
 E: enquiries@drivemedical.co.uk
 W: www.drivemedical.co.uk/technical

TITLE:
**POWERSTROLL (FROM 2010)
 S-DRIVE PARTS MANUAL**

DATE: 08-09-10
 DRAWN: PMK
 SCALE: NTS

SHEET No: 1 of 1
 REVISION No: C
 DWG No: Z24098