



Code	Description	Specification
Z16215-24	TR39 Seat Canvas	Nylon, 47cm (W) x 41cm (D), 4x2 holes at 5, 15,
		25, 35cm
Z16265-24	TR39 Back Canvas	Nylon, 47cm (W) x 48cm (D)
Z12325-24SV	TR39 Push Handle	Silver Vein
Z12350-24	TR39 LHS Half Fold Mechanism	Silver Vein
Z12360-24	TR39 RHS Half Fold Mechanism	Silver Vein
Z12410-24	TR39 Front Castor Wheel & Fork	Black Fork, shaft 12.7mm (1/2") dia x 60mm
		8-spoke wheel, 20cm x 2.8cm black tyre & hub,
		Wheel bearings for 8mm axle.
Z16460-24	TR39 Rear Wheel	8-spoke wheel, 20cm x 2.8cm black tyre & hub,
		Wheel bearings for 11mm axle.
Z11592-24	TR39 Arm Pad	Full length vinyl, 25.5cm hole spacing.
		35.5cm x 5.6.
Z12600-24SV	TR39 RHS Leg Rest Ass'y	Silver Vein with black inner stem, 38mm pin
		height spacing
Z12650-24SV	TR39 LHS Leg Rest Assy	Silver Vein with black inner stem, 38mm pin
		height spacing
Z12700-24	TR39 RHS Brake Unit	Black pressed steel, mtd in front of back post
Z12750-24	TR39 LHS Brake Unit	Black pressed steel, mtd in front of back post
Z1802-24	TR39 Rubber Stepper Pad	6.5cm (L) x 3.3cm (W), to fit 22mm dia tubing
Z11024-24	TR39 Plastic Seat Guide (Front)	See guide
Z11025-24	TR39 Plastic Seat Guide (Rear)	See guide
Z15900-01A	Lap Belt	2-piece lap belt for aluminium wheelchairs
Z12100-24	Hub Cap	Rear wheel hub cap TR-39ESV

Document Issued 09/15