

Parts List for Lift Chairs – Instore and Private Label Chairs

(Issued Feb-19)

1. Single Motor Lift Chairs – Fabric and PU (CLR25, CLR30 and CLR31)

Parts available below (chairs produced upto Jan 2019)

Code	Description	Spec
C10-203-00100	Power Supply	<i>Input: 100-240V 50/60Hz 1.6A UK plug figure-8 lead Output: 29V 1.8A 2-pin speaker socket</i>
C10-203-00200	Handset	<i>2-button handset with coiled cable and 5-pin 180° DIN socket</i>
C09-203-00100	Actuator	<i>Richmat HJA58 - Stroke 173mm / Install 397mm DC 24V IP20 with 5-pin 180° DIN plug and 2-pin speaker plug</i>

Parts available below (chairs produced after Feb 2019)

Code	Description
C09-203-00400	Actuator for Single motor - CLR 25-30-31 (Wi-Base 2019 onwards)
C10-205-00200	Handset, 2-button (Wi-Base 2019 onwards)
C10-205-00300	PSU w/o battery backup (Wi-Base 2019 onwards)

2. Dual Motor Lift Chairs – Fabric and PU (CLR26, CLR32, CLR33 and CLR38)

Parts available below (chairs produced up to Jan 2019):

Code	Description	Spec
C10-203-00100	Power Supply	<i>Input: 100-240V 50/60Hz 1.6A UK plug figure-8 lead Output: 29V 1.8A 2-pin speaker socket</i>
C10-203-00400	Handset	<i>4-button handset with coiled cable and 5-pin 180° DIN socket</i>
C09-203-00200	Actuator	<i>For legrest</i>
C09-203-00300	Actuator	<i>For backrest</i>

Parts available below (chairs produced after Feb 2019)

Code	Description
C09-203-00400	Actuator for Single motor - CLR 25-30-31 (Wi-Base 2019 onwards)
C09-203-00500	Actuator for Dual motor(Riser) - CLR 26-32-33 (Wi-Base 2019 onwards)
C09-203-00600	Actuator for Dual motor(Recliner) - CLR 26-32-33 (Wi-Base 2019 onwards)
C10-205-00100	Handset, 4-button (Wi-Base 2019 onwards)
C10-205-00200	Handset, 2-button (Wi-Base 2019 onwards)
C10-205-00300	PSU w/o battery backup (Wi-Base 2019 onwards)